Dispositif médical de classe IIa

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Medical device Class IIa

Indications: removal of the superficial layers of the skin in the treatment of the epidermal affections of the foot

Contraindications: allergy to stainless steel

Instruction for use:

1. Check the expiry date indicated on the box and on the individual packaging
2. On first use, cut along the dotted lines the location of the dispenser
3. Take an individual package containing a sterile blade
4. Check the integrity of the primary packaging
5. Open the primary packaging, without touching the cutting part, to release the attachment part
6. Fix the blade to a suitable handle
7. Proceed to the patient care immediately after blade installation
8. After care, remove the blade from the handle with nail clippers
9. Dispose immediately after use according to the national provisions in force

WARNINGS

Sharp medical device intended for professionals trained in its use.

Store in a dry and ventilated place. Do not expose to direct sunlight. Do not wet. Do not crush the packaging. Store away from corrosive substances.

The blades must be kept in their original box and in their individual packaging until use.

In the following cases, do not use the blade:

- the expiry date has passed,
- the packaging is damaged or unintentionally opened before use

Single-use medical device, sterilized by Gamma rays. Do not re-sterilize. Risk of infection if reused.

To avoid any risk of injury, the blade must be attached to a handle before use (EN 27740).

Any serious incident occurring in connection with the device must be notified immediately to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established

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COBLENTZ MEDICAL BLADES INDUSTRY (CZ.MB.I) – Le Courbet – 531 route du Bellay 49650 ALLONNES - S.A.S AU